

EXHAUST ARM



Articulated 80Φmm Arm

ADCWAH-80

Arms are movable and are excellent point of use exhaust control.



Transparent (PC)

Non-Transparent (PP)

Model	ADCWAH-80-15	ADCWAH-80-18	ADCWAH-80-21	ADCWAH-80-24
Dimensions	80φ x 1500mm	80φ x 1800mm	80φ x 2100mm	80φ x 2400mm
Material	Polypropylene (PP)			
Bracket Material	Industrial grade electro galvanised steel that is powder coated			



80Φ mm diameter



Damper control



360 degree rotation



Built in carbon filter



Easy disassembly, assembly and moving with the use of control bolt



Quality high tension spring

Dome Selection

- 400Φmm Transparent PC Dome (Standard)
- 400Φmm Non-transparent PP Dome
- 200Φmm PP Dome



Articulated 110Φmm Arm

ADCWAH-110

Arms are movable and are excellent point of use exhaust control.



Transparent (PC)

Non-Transparent (PP)

Model	ADCWAH-110-18	ADCWAH-110-21	ADCWAH-110-24
Dimensions	110φ x 1800mm	110φ x 2100mm	110φ x 2400mm
Material	Polypropylene (PP)		
Bracket Material	Industrial grade electro galvanised steel that is powder coated		



110Φ mm diameter



Wind volume controllable damper



Easy handling by user with
360 degree rotation



Better with 450Φ dome



Easy disassembly, assembly and
moving with the use of control bolt



Quality gas spring

Dome Selection

- 450Φ Transparent PC Dome (Standard)
- 450Φ Non-Transparent PP Dome



Flexible Hose Arm Hood

ADCWMH-60

Flexible Hose, can be connected to any exhaust system

Model	ADCWMH-60
Dimensions	77Φ x 60mm
Outer Diameter	77Φ mm
Inner Diameter	75Φ mm
Material	POM



Designed to minimize air stagnation



Easy adjustment of length



Working dynamics



Resistant against chemical, heat, weather and abrasion



Dome Selection

- 250Φ Dome
- 150Φ Dome

250Φ Dome

150Φ Dome



Singapore

Advancelab (S) Pte Ltd

253 Kaki Bukit Ave 1,
Singapore 416061
Tel: +65 6448 8255
Fax: +65 6448 9833
Email: info@advancelab.com.sg

Malaysia

Advancelab Sdn Bhd

No. 3388, Jalan Pekeliling Tanjung
27/2, Kawasan Perindustrian Indahpura,
81000 Kulaijaya, Johor, Malaysia.
Tel: +607 660 8877
Fax: +607 660 8866
Email: info-my@advancelab-global.com

Myanmar

Advancelab Scientific &
Engineering Co., Ltd

No.(81/1)-1A, Myin Thar 7th Street,
(14/1)Ward, South Okkalapa Township,
Yangon, Myanmar.
Tel: +95 (9) 779753802
Email: info-mm@advancelab-global.com

Indonesia

PT. Advancelab Saintifik

253 Kaki Bukit Ave 1,
Singapore 416061
Tel: +65 6448 8255
Fax: +65 6448 9833
Email: info@advancelab.com.sg

Thailand

Prima Scientific
Co., Ltd.

147 / 170 - 171 Baromrajchonnane
Road, Arunamarin Bangkoknoi, Bangkok
10700, Thailand.
Tel: +66 2 884 9480
Fax: +66 2 884 6441
Email: primasci@primasci.com

U.A.E (Dubai)

Advancelab FZCO

253 Kaki Bukit Ave 1,
Singapore 416061
Tel: +65 6448 8255
Fax: +65 6448 9833
Email: info@advancelab.com.sg

Vietnam

M&T International Trading
Services Co., Ltd.

46A Phan Dinh Phung, Tan Thanh Ward,
Tan Phu District, Ho Chi Minh City
700000, Vietnam.
Tel: +84 28 6656 0610
Email: info@advancelab.com.sg

