

# NEW JERSEY



BY ADVANCELAB







**BUILT TOUGH**





100  
876

450



1190  
876

914


457

914 / 1473


- 01 Reagent Rack
- 02 Table Top
- 03 Drawer
- 04 Cabinet
- 05 Knee Space
- 06 Adjustable Shelves
- 07 Service Dropper

# Cut Sheet


## Wall Cabinet




## Mobile Pedestal



## Sink Unit



## Corner Cabinet





# COLOUR

- 01 Ocean Blue
- 02 Light Blue
- 03 Pastel Green
- 04 Emerald Green
- 05 Yellow
- 06 Yellow Orange
- 07 Traffic Red
- 08 Purple White
- 09 White
- 10 Cream
- 11 Light Grey
- 12 Dark Grey
- 00 Black

# Accessories

## Taps



**1 Way Water Faucet**  
 ADBM104 - 3B1  
 Bench Mount, Swivel Neck & Fixed Nozzle.  
 Weight: 1.18 kg



**2 Way Water Faucet**  
 ADBM110 - 3B1  
 Bench Mount & Fixed Nozzle.  
 Weight: 1.45 kg



**3 Way Water Faucet**  
 ADBM106 - 5B1  
 Bench Mount, Swivel Neck & Fixed Nozzle  
 Weight: 2.18 kg



**1 Way Bench Gas Valve**  
 ADOM010 NG  
 Lift & Turn Operation & Fixed Nozzle  
 Weight: 0.65 kg




**2 Way 90° Bench Gas Valve**  
 ADOM020 NG  
 Lift & Turn Operation & Fixed Nozzle  
 Weight: 0.83 kg



**2 Way 180° Bench Gas Valve**  
 ADOM025 NG  
 Lift & Turn Operation & Fixed Nozzle  
 Weight: 0.81 kg

## Reagent Rack



Reagent Rack

## Electric Sockets



**Single 3-Pin Power Socket**  
 Electric Power Socket  
 13A



**Twins 3-Pin Power Socket**  
 Electric Power Socket  
 13A

## Drying Rack



**AD 450 Drying Rack**  
 High impact Polyethylene drying rack. The back plate of the 72 place support is made from a single mould of 4 mm rigid HIPS.  
 Size: 450mm x 630mm x 110mm

## Safety Equipment



**Single Bench Mount Eyewash**  
 ADE200SS  
 Single Nozzle Bench Spray Unit, Vertical Extraction (Top Mounting)



**Twins Bench Mount Eyewash**  
 ADE300SS  
 Twins Nozzle Bench Spray Unit, Vertical Extraction (Top Mounting)

## Dilution Tank



**AD250 PPDT**  
 Advancelab dilution tank made by full polypropylene construction. Available in 250mm diameter.  
 Size: 250mm x 300mm x 250mm

## Worktops

	Type	Thickness	Colour
	Phenolic Resin	16mm (Standard) 25mm	Grey
	Epoxy Resin	19mm (Standard) 25mm	Grey, Black
	Stainless Steel 304 Stainless Steel 316	1.2mm (Standard) 1.5mm	





## Singapore

### Advancelab (S) Pte Ltd

253 Kaki Bukit Ave 1, Singapore 416061  
Tel: +65 6448 8255  
Fax: +65 6448 9833  
Email: info@advancelab.com.sg  
Website: www.advancelabglobal.com

## Malaysia

### Advancelab Sdn Bhd

No. 3388, Jalan Pekeliling Tanjung 27/2,  
Kawasan Perindustrian Indahpura,  
81000 Kulajaya, Johor, Malaysia.  
Tel: +607 660 8877  
Fax: +607 660 8866  
Email: info-my@advancelab.com.sg

## Myanmar

### Advancelab Scientific & Engineering Co., Ltd

No.(81/1)-1A, Myin Thar 7th Street,  
(14/1)Ward, South Okkalapa Township,  
Yangon, Myanmar  
Tel: +95 (9) 779753802  
Email: info-mm@advancelab-global.com

## Indonesia

### PT Advancelab Saintifik

253 Kaki Bukit Ave 1, Singapore 416061  
Tel: +65 6448 8255  
Fax: +65 6448 9833  
Email: itnfo@advancelab.com.sg

## Thailand

### Prima Scientific Co., Ltd.

147 / 170 - 171 Baromrajchonnanee Road,  
Arunamarin Bangkoknoi, Bangkok 10700,  
Thailand.  
Tel: +66 2 884 9480  
Fax: +66 2 884 6441  
Email: primasci@primasci.com

## U.A.E (Dubai)

### Advancelab FZCO

No.FZJOAB1310 Jafza One, Jebel Ali  
Dubai, United Arab Emirates  
Email: info@advancelab.com.sg

NEW  
JERSE