

PASSBOX

BY ADVANCELAB






Advancelab offers stainless steel (316 & 304) as well as powder coated electrogalvanised steel passboxes.

- Double-layered design
- Mechanical interlock is reliable, maintenance-free, fail safe and maintains cleanroom integrity by preventing both doors from being opened at same time
- Tempered glass provides a clear view of the internal chamber
- Designed simply for years of hassle-free use


Cut Sheet

Mechanical Interlock



- 01 Hinges
- 02 Glazed Glass Window
- 03 Handle with Lock
- 04 Interlocking System
- 05 Wall Frame (optional)

Wall Frame (optiona)



Model	AD01PB	AD02PB	AD03PB
External Dimensions (W x D x H)	400 x 460 x 600mm	600 x 460 x 600mm	900 x 460 x 900mm
Internal Work Area Dimensions (W x D x H)	300 x 400 x 475mm	500 x 400 x 475mm	800 x 400 x 775mm
Material	Electrogalvanised Steel / Stainless Steel 304 / Stainless Steel 316		
Window	Glazed Glass / No Window		
Opening	Left Hand Open / Right Hand Open		
Wall Frame (optional)	Electrogalvanised Steel / Stainless Steel 304 / Stainless Steel 316		
Model with Wall Frame	AD01PB-WF	AD02PB-WF	AD03PB-WF

Singapore

Advancelab (S) Pte Ltd

253 Kaki Bukit Ave 1
Singapore 416061
Tel: +65 6448 8255
Fax: +65 6448 9833
Email: info@advancelab.com.sg

Malaysia

Advancelab Sdn Bhd

No. 3388, Jalan Pekeliling Tanjung 27/2,
Kawasan Perindustrian Indahpura,
81000 Kulaijaya, Johor, Malaysia.
Tel: +607 660 8877
Fax: +607 660 8866
Email: info-my@advancelab-global.com

Myanmar

Advancelab Scientific &
Engineering Co., Ltd

No.(81/1)-1A, Myin Thar 7th Street,
(14/1)Ward, South Okkalapa Township,
Yangon, Myanmar
Tel: +95 (9) 779753802
Email: info-mm@advancelab-global.com

U.A.E (Dubai)

Advancelab FZCO

253 Kaki Bukit Ave 1
Singapore 416061
Tel: +65 6448 8255
Fax: +65 6448 9833
Email: info@advancelab.com.sg

Indonesia

PT Advancelab Saintifik

253 Kaki Bukit Ave 1
Singapore 416061
Tel: +65 6448 8255
Fax: +65 6448 9833
Email: info@advancelab.com.sg

Thailand

Prima Scientific Co., Ltd.

147 / 170 - 171 Baromrajchonnane
Road, Arunamarin Bangkoknoi, Bangkok
10700, Thailand.
Tel: +66 2 884 9480
Fax: +66 2 884 6441
Email: primasci@primasci.com

