

# AdvancelabPro® Steel

CASEWORK



# Our Story

**Island**

**Bench**

**Ergonomics**

**Construction**

**Cut Sheet**

**Palette**

**Accessories**

---

**Taps**

**Reagent Rack**

**Electric Sockets**

**Safety Equipment**

**Dilution Tank**

**Worktop**



A South-East Asian and beyond **MANUFACTURER** with HQ in Singapore and manufacturing base in Malaysia.

Our DNA consists of forward thinking designs, **FLEXIBLE** layouts and creativity.

We are **EXPORTING** to all continents and count more than 40 regular distributors who have been with us for years.

We build the correct **RELATIONSHIPS**.

Our large range of manufacture include everything from pass boxes to cleanbooths, allowing us to give **SOLUTIONS** across the board.





# PEOPLE COMPANY PRODUCTION COUNTRY

from  
SINGAPORE  
to the  
world




## Island



## Bench



Ergonomics






## Construction

**FLEXIBILITY** - standard modules of 450 and 900 mm allow us to set up a lab and easily reconfigure it to suit future needs.

**LONG LASTING** electrogalvanised steel construction with 1.0 mm thickness with a heavy powder coated finish for a lush look that will last years in robust laboratory environments.


**SLIDING KNEE SPACES** - this makes it incredibly flexible on site, it can be adjusted at site to suit local construction conditions. A very forgiving set up.






Cut Sheet


Underbench




Sink Cabinet



Corner Cabinet



Wall Hung



# Palette



Light Blue  
PR12/363/CS21

Ocean Blue  
PR12/3106/CS6

Emerald Green  
PR12/2045/CS2  
Pastel Green  
PR12/227/CS13

Signal Yellow  
FF160/164/C

Orange  
FF160/6178/C4

Traffic Red  
PR12/526/CS28

Black  
PR12/708/CS

Dark Grey  
PR12/892/CS10

Light Grey  
PR12/808/CS10

Cream  
PR12/4249/CS1

White  
LB1000

Purple White  
FF160/12/C1



# Accessories

## Taps



**1 Way Water Faucet**  
ADBM104 - 3B1  
Bench Mount, Swivel Neck & Fixed Nozzle.  
Weight: 1.18 kg



**2 Way Water Faucet**  
ADBM110 - 3B1  
Bench Mount & Fixed Nozzle.  
Weight: 1.45 kg



**3 Way Water Faucet**  
ADBM106 - 5B1  
Bench Mount, Swivel Neck & Fixed Nozzle  
Weight: 2.18 kg

## Electric Sockets



**Single 3-Pin Power Socket**  
Electric Power  
Socket 13A



**Twins 3-Pin Power Socket**  
Electric Power  
Socket 13A



**AD 450 Drying Rack**  
High impact Polyethylene drying rack. The back plate of the 72 place support is made from a single mould of 4mm rigid HIPS.  
Size: 450mm x 630mm x 110mm

## Safety Equipment



**1 Way Bench Gas Valve**  
ADOM010 NG  
Lift & Turn Operation & Fixed Nozzle  
Weight: 0.65 kg



**2 Way 90° Bench Gas Valve**  
ADOM020 NG  
Lift & Turn Operation & Fixed Nozzle  
Weight: 0.83 kg



**2 Way 180° Bench Gas Valve**  
ADOM025 NG  
Lift & Turn Operation & Fixed Nozzle  
Weight: 0.81 kg




**Single Bench Mount Eyewash**  
ADE200SS  
Single Nozzle Bench Spray Unit, Vertical Extraction (Top Mounting)



**Twins Bench Mount Eyewash**  
ADE300SS  
Twins Nozzle Bench Spray Unit, Vertical Extraction (Top Mounting)

## Reagent Rack



## Worktops

	Type	Thickness	Colour
	Phenolic Resin	16mm (Standard) 25mm	Grey
	Epoxy Resin	19mm (Standard) 25mm	Grey, Black
	Stainless Steel 304 Stainless Steel 316	1.2mm (Standard) 1.5mm	

## Dilution Tank



**AD250 PPDT**  
Advancelab dilution tank made by full polypropylene construction. Available in 250mm diameter.  
Size: 250mm x 300mm x 250mm







## Singapore

### Advancelab (S) Pte Ltd

253 Kaki Bukit Ave 1, Singapore 416061  
Tel: +65 6448 8255  
Fax: +65 6448 9833  
Email: info@advancelab.com.sg

## Malaysia

### Advancelab Sdn Bhd

No. 3388, Jalan Pekeliling Tanjung 27/2,  
Kawasan Perindustrian Indahpura,  
81000 Kulaijaya, Johor, Malaysia.  
Tel: +607 660 8877  
Fax: +607 660 8866  
Email: info-my@advancelab-global.com

## Myanmar

### Advancelab Scientific & Engineering Co., Ltd

No.(81/1)-1A, Myin Thar 7th Street,  
(14/1)Ward, South Okkalapa Township,  
Yangon, Myanmar.  
Tel: +95 (1) 572393 / +95 (9) 779753802  
Fax: +95 (1) 572393  
Email: info-mm@advancelab-global.com

## Indonesia

### PT Advancelab Saintifik

253 Kaki Bukit Ave 1, Singapore 416061  
Tel: +65 6448 8255  
Fax: +65 6448 9833  
Email: info@advancelab.com.sg

## Vietnam

### M&T International Trading Services Co., Ltd

46A Phan Dinh Phung, Tan Thanh Ward,  
Tan Phu District, Ho Chi Minh City  
700000, Vietnam.  
Tel: +84 28 6656 0610  
Email: info@advancelab.com.sg

## Thailand

### Prima Scientific Co., Ltd.

147 / 170 - 171 Baromrajchonnane Road,  
Arunamarin Bangkoknoi, Bangkok 10700,  
Thailand.  
Tel: +66 2 884 9480  
Fax: +66 2 884 6441  
Email: primasci@primasci.com

## U.A.E (Dubai)

### Advancelab FZCO

253 Kaki Bukit Ave 1, Singapore 416061  
Tel: +65 6448 8255  
Fax: +65 6448 9833  
Email: info@advancelab.com.sg